Mittelmann

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
RsqTec **COCOON**

DIN EN 1497:2007

DIN EN 1865-1:2015 Ziffer 5.4.1

RSQTEC
COCOON

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DE

Instructions for use / test book

GB

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Kennnummer: CE 0158

Instructions for Use

These instructions for use are general instructions for the product described and are not a substitute for user training. Intensive instruction and practice in the use and maintenance of the rescue stretcher with carrying function is mandatory, in addition, the user must also be both physically and mentally fit. Inadequate training, incorrect use or misuse of the product can lead to accidents. The advice and instructions contained in the instructions for use must be carefully taken note of and complied with.

Foreword

The strengths of the Development department and field service specialists of Mittelmann Sicherheitstechnik have been brought together for the development of the "COCOON" rescue stretcher with carrying function. Hence, it was possible to combine development know-how for special rescue equipment together with the industry-specific requirements and ideas from real-life application.

During the development period, the rescue stretcher with carrying function was tested in a wide variety of application scenarios during ongoing training operations and was continuously further developed.

1. General Application

The rescue stretcher with carrying function has been developed for the lying rescue transport of injured persons. The rescue system serves for the horizontal and vertical transport of injured people from exposed locations for rope and, especially, air rescue. The enclosing of the person over the entire length takes place by means of tear-resistant and water-repellent nylon fabric and the injured person is secured with an integrated rescue harness. The objective is to ensure the safety of the person and to minimize the physical exertion of the rescue workers as much as possible. The rescue stretcher with carrying function must only be used within the specified operating conditions and for its intended purpose.

Through interaction, every combination of different equipment can give rise to unforeseen hazardous situations and negatively affect the safety of the user.

2. Application Areas

Mountain, aerial and high-altitude rescue, for example, from areas that are difficult to access or from impassable terrain.

3. User Training

The emergency personnel must be adequately trained and be familiar with the use of the rescue stretcher with carrying function before its first use. In addition, in the introduction and the subsequent further training, the repeated teaching of these instructions for use is necessary.

The training course must be strictly documented and be repeated at least once per year. The type, the scope and the date of the training must be chronologically recorded.

4. Product Description


The "COCOON" rescue stretcher with carrying function is a combination of a rescue harness and a patient transport device and has been tested and certified in accordance with DIN EN 1497:2007 and DIN EN 1865-1:2015 Paragraph 5.4.1.

The system is made of textile materials for the rescue and transport of the injured person. The textile sheath (a) encases the entire length of the person. Fastening is done by means of large size Velcro strips. In the upper part of the RsqTec COCOON, there is a sight-free head cover (b), opposite this is the foot cover (c). Four lashing belts (d) with adjustable snap buckles are used to secure the protective tray close to the body of the person. A system of supporting cross and diagonal belts that bear the body weight of the person during the rescue and transport runs inside the RsqTec COCOON. The inner harness consists of an ordinary rescue harness that is constructed as follows: The shoulder straps (e) are length-adjustable and are connected to the horizontal chest strap (f) by means of two snap buckles. Additional supporting straps are the leg straps right (g) and left (h) that likewise can be adjusted for each person by means of the snap buckles. All five pairs of the carrying straps (i) have been designed for the lifting or transferring of the RsqTec COCOON in an inclined or vertical position but only one pair is used at a time. The choice of the carrying strap pairs in conjunction with the suspension loop (k) and two snap hooks is dependent on the respective use and rescue case. The rescue and / or transport of the person lying horizontally is carried out via the 10-point suspension loops (j) that are equipped with two triangular locking rings. The rescue of the person suspended vertically takes place with the aid of the head loops green and red (l). In addition, there are carrying straps with handle functions (i) on the outside of the RsqTec COCOON for manual transport.

5. Designation of the Individual Parts

No.	Designation	Attachment Points for the Following Lying Positions
(a)	Textile sheath (red)	
(b)	Sight-free head cover (red)	
(c)	Foot cover	
(d)	Lashing belts (black) 4x	
(e)	Shoulder straps (blue)	
(f)	Chest strap (blue)	
(g)	Leg strap right (green)	
(h)	Leg strap left (red)	
(i)	Carrying strap pairs ((1-5), 1x yellow, 4x orange)	Vertical and inclined lifting (ONLY in combination with (k))
(j)	Delta lock ring with black shackle 2x (10-point suspension (black))	Horizontal transport Horizontal lifting
(k)	Suspension loop (yellow)	Vertical transport Vertical lifting
(l)	Head loops green and red	Vertical transport (ONLY in combination with (k))

5.1 Illustrative Representation


6. Preparation and Belt Fitting of the RsqTec COCOON


The preparation of the RsqTec COCOON is carried out as follows:

The rolled up and packed RsqTec COCOON is pulled out of the transport bag. Next, the fastening strap with the clamp lock is released so that the RsqTec COCOON can be unfolded lengthwise.



After the RsqTec COCOON is fully unfolded, all four corners of the protective tray are passed through the respective slotted holes in the order shown in the picture sections and connected and tightened via the snap buckles.

The preparation of the RsqTec COCOON is completed by inserting the two sail battens into the pockets provided. The pockets are located on the inside of the two textile sheaths on the sides.

Before starting the belt fitting, all belts and snap buckles are opened and placed to the side. Now the injured person is placed with his / her back on the RsqTec COCOON so that the top of the head is approx. 5 cm below the inner edge of the protective tray. Then the step-by-step instructions for correct belt fitting, presented in table form, are carried out.

<i>Illustration</i>	<i>Description</i>
	Pass the long side of the leg straps (red (h) and green (g) from Section 5.1) from the side under the respective thigh of the person. Lead up into the crotch and connect with the corresponding snap buckles (the snap buckles must engage audibly and tangibly). The strap is tightened (not constricting) without changing the position of the person.
	The head fixation (from year of manufacturing 01.01.2023) consists of an adhesive textile fabric and is fixed to the forehead of the injured person or to the front surface of the helmet with the red and green webbing. The position of the head fixation can be adjusted using the flexible black webbing.

Place the shoulder strap ((e) from chapter 5.1) lengthwise over the head of the person, guide it around the shoulders and connect it in front of the chest with the side parts of the chest strap ((f) from Chapter 5.1). The strap is tightened (**not constricting**). It is imperative that the full freedom of movement of the chest for breathing is maintained.

Connect the protective tray with the aid of the four lashing belts ((d) from Chapter 5.1) and the snap buckles and tighten them tightly around the body of the injured person (**not constricting**).

Place the textile sheath ((a) from Chapter 5.1) snugly around the person and close it with the help of the Velcro fastener.

The next step is to close the sight-free head cover ((b) from Chapter 5.1) by connecting three press studs on the plastic protective tray. The foot cover ((c) from Chapter 5.1) is attached in the same way, in doing so, it must be ensured that the zip fastener is closed as shown in the illustration.

The injured person is now prepared for horizontal and vertical transport with the RsqTec COCOON.

Optionally Available Accessories	Description
	Foot sling for RsqTec TURTLE and COCOON for leg relief during vertical personal rescue, reduces the risk of orthostatic shock. The attachment of the foot sling is made by means of snap buckles on the black webbing loops below the leg straps.
	The face shield covers the entire face with a flexible visor and protects the injured person from various hazards. These can be of a mechanical, thermal, chemical or optical nature. Adaptively attachable to the RsqTec COCOON from year of manufacture 01.01.2023.
	The anti-rotation system is operated via the ergonomic handle and allows you to control or reduce the rotation effects that occur during recovery with the helicopter (downwash effect). Adaptively attachable to the RsqTec COCOON from year of manufacture 01.01.2023.
	The neck support pillow is intended to relieve tension and is placed on the black fleece band in the head area of the rescue stretcher and thus connected. Adaptively attachable to the RsqTec COCOON from year of manufacture 01.01.2023.

	Waterproof transport bag, Ø 460x1100mm, with carrying handle and 2 shoulder straps (rucksack with clamp locks), blue/black with transparent window and Mittelmann logo. Storage pocket for DIN A5, extra loops.
	Carrying mattress: Multi-chamber vacuum mattress 2000x800 mm with 6 carrying handles (4.7 kg). Including plastic pump.

7. Application Possibilities under Different Situation Positions

The following possible applications only take place after the RsqTec COCOON has been successfully set up and the belts have been correctly fitted to the person (see Chapters 6 and 7).

The attachment point for all positions must lie above the user and guarantee a minimum load carrying capacity of 10 kN.

Horizontal lifting and transport of the person with the RsqTec COCOON.

- (i) 2x Delta lock ring with black shackle (10-point suspension).

The two Delta lock rings (10-point suspension) ((j) from Chapter 5.1) are taken out of their side pannier and attached to the connection device of the lifting gear.

The injured person is lifted with the aid of a suitable hoist, whereby the horizontal position of the person must be constantly checked and, if necessary, be corrected by changing the lying position of the person.

Vertical or inclined lifting of the person with the RsqTec COCOON

Note: In order to prevent injury to the person, the person must only be lifted in the vertical position for a short time (for a few minutes).

- (i) Connect the suspension loop (k) with one carrying strap pair (i).

The connection device of the lifting gear is attached to the suspension loop (k) which is connected with two snap hooks to a carrying strap pair ((i) from Chapter 5.1). Which of the carrying strap pairs is selected in conjunction with the suspension loop (k) is dependent on the respective application and rescue case.

During this process, special attention must be paid so that no excessive pressures occur for the person. Where necessary, the position of the relevant belt must be checked.

Vertical transport of the person with the RsqTec COCOON

Note: In order to prevent injury to the person, the person must only be lifted in the vertical position for a short time (for a few minutes).

- (i) Connect the suspension loop (k) with the head loops green and red (l).

The connection device of the lifting gear is attached to the suspension loop (k) which is connected with two snap hooks to the head loops green and red ((l) from Chapter 5.1). In order to be able to transport the injured person in a vertical position with the aid of a suitable lifting device, the green and red head loops must be located outside the protective tray and passed through the openings provided for this purpose.

Recommendation:

The additional article ("RsqTec Foot Sling TURTLE/COCOON" – Article Code: 200145) is for RsqTec TURTLE and COCOON leg relief during vertical personal rescue, reduces the risk of orthostatic shock.

The foot sling is attached to the black webbing loops below the leg straps using snap buckles (see the following picture instructions).

8. Technical Data and Parameters

Physical Dimensions

Maximum length	2200 mm = 2.2 m
Maximum width	550 mm = 0.55 m
Maximum height	340 mm = 0.34 m
Diameter (packed size)	Ø = 420 mm = 0.42 m
Height (packed size)	880 mm = 0.88 m
Weight without person	ca. 11.5 kg

Load Limits


Maximum nominal load	140 kg
----------------------	--------

Temperature Limits

Minimum outdoor temperature	- 40°C = - 40°F
Maximum outdoor temperature	+ 60°C = + 140°F

9. Meaning of the Marking

Mittelmann Sicherheitstechnik	Manufacturer Mittelmann Sicherheitstechnik GmbH & Co. KG
Bessemerstraße 25, 42551 Velbert	Postal address of the manufacturer
www.mittelmann.com	Internet address of the manufacturer
RsqTec	Product designation
COCOON	Type designation
DIN EN 1497:2007 DIN EN 1865-1:2015 Paragraph 5.4.1	Number and year of the documents with which the equipment complies
Nominal load: 140 kg	Specification of the maximum permissible body weight
Made in Germany	Designation of origin
CE 0158	Identification number of the quality monitoring notified inspection body
Year of manufacture: xx/xxxx	Month / year of manufacture
	Symbol for indicating that the instructions for use must be observed
Serial No.: xxxxxxxxxx	Serial number

10. Material Specifications

RsqTec *COCOON* consisting of the following materials:

- Straps in polyamide.
- Material sections in polyamide.
- Sewing thread in polyester.
- Fitting components in aluminium, steel and stainless steel.
- Plastic components in polyethylene.
- Snap buckles in aluminium and stainless steel.
- Plastic parts in plastic material.
- Carabiner hooks in aluminium or stainless steel.

11. Storage and Transport

- The RsqTec *COCOON* must be stored dry, protected from direct sunlight and at room temperature (approx. 17 °C – 23 °C).

- At the storage location, all types of contact with acids, corrosive liquids and oils must be prevented.

- Whenever the RsqTec *COCOON* becomes wet, it must be dried naturally before being stored.
- During the transport of the RsqTec *COCOON*, it must be stored in a suitable transport bag.

12. Damage Analysis

Before each use of the rescue stretcher, it must be checked for the following signs of wear and / or damage:

- Strap or textile damage due to abrasion, cuts, heat or chemicals.
- Damage to any type of seam.
- All types of damage to the load-bearing textile structure.
- All types of contamination through technical liquids that can no longer be cleaned.
- Damage to the Velcro fasteners.
- Deformation or corrosion of the fitting components or snap buckles.
- Defective snap buckles or clasps
- Check the legibility of the product marking.

Immediately one of the damage analysis points described above is detected before use, the rescue stretcher with carrying function must be taken out of use.

13. Safety Requirements

- Through interaction, every combination of different equipment can give rise to unforeseen hazardous situations and negatively affect the safety of the user.
- All instructions for the protection of the user and the rescue stretcher with carrying function given in these instructions for use must be strictly observed. The rating plate must always be visible and legible.
- In order to avoid a suspension trauma, an emergency plan on the part of the user must be available for every type of use of the rescue stretcher with carrying function. A suspension trauma is a life-threatening state of shock that occurs after long, motionless hanging in a harness system. This can lead to unconsciousness and to cardiac arrest.
- This emergency plan must ensure that the vertical hanging time of the person in the rescue stretcher is kept as short as possible and so prevent massive damage to health.
- It is vital that the belt fitting of the rescue stretcher to the person is regularly checked during its use.

- For the safety of the person to be rescued, it is essential that before every use of the rescue stretcher with carrying function the free space necessary for the rescue operation is ensured so that no impact on the ground or with other obstacles is possible.
- Hazards, such as for example, extreme temperature, stress from sharp edges, cuts, abrasion, slack rope formation, chemical effects, electrical influences or climatic effects can adversely affect the function of the equipment.

14. Inspection and General Instructions

- Should there be doubts about the safe condition of the system or the component parts, these must be replaced immediately. This must be carried out by the manufacturer or another qualified person.
- Systems stressed by a fall must be taken out of use and be returned to the manufacturer or a qualified workshop for servicing and re-examination.

- **Do not make any changes to the product.**

- **Repairs must only be carried out by the manufacturer.**

- The instructions for use / test book are delivered with every rescue stretcher with carrying function. These must be kept with the equipment.


- Strictly observe the instructions for use of the products used in connection with the rescue stretcher with carrying function as well as the relevant applicable accident prevention regulations.
- The Declaration of Conformity can be downloaded from the www.mittelmann.com page.
- As required, but at least once per year, an inspection by the manufacturer or a qualified person must be carried out (see also Point 12 Damage Analysis). This inspection with the date must be documented in the test book supplied. In addition, the next inspection date must be clearly recorded on the rating plate.
- The rating plate with the year of manufacture and the serial number can be found in the rating plate pocket on the left shoulder strap (e).
- If the equipment is sold on in another country, the re-seller must provide instructions for the use, the maintenance, the regular inspection and repair in the language of the other country.
- The rescue stretcher with carrying function must not be used for fall arrest purposes.

15. Cleaning Instructions

- The rescue stretcher with carrying function must be thoroughly cleaned after every use and be stored dry in an airy and shady space.

- Avoid contact with acids, corrosive liquids and oils. If unavoidable, clean as soon as possible by hand washing with mild detergent and plenty of water (approx. 30°C). Always allow to dry naturally, never in the vicinity of fire or other sources of heat.

- Please do not use dry cleaning, tumble dryers, spin dryers or bleaching processes.

- For questions regarding disinfection, please contact the manufacturer.
- Any body fluids on the rescue stretcher with carrying function must be cleaned off using water only.

16. Service Life

The service life is dependent on the respective usage conditions. On the basis of past experience, under normal usage conditions the service life of the **COCOON** rescue stretcher with carrying function can be assumed to be 10 years. The year of manufacture of the rescue stretcher is recorded on the rating plate.

In the event of any kind of damage to any component, the rescue stretcher with carrying function must be available to the manufacturer for the evaluation of the damage. The continued use of the rescue stretcher with damaged components shall lead to the immediate loss of any and all rights to liability claims.

Notified body involved in the production monitoring:

DEKRA Testing and Certification GmbH

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