Mittelmann

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RsqTec TURTLE FrameX

DIN EN 1497:2007
DIN EN 1865-1:2015

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DE

Instructions for use/ test book

GB

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FR

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DEKRA Testing and Certification GmbH

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Instructions for use

These instructions for use are general instructions for the product described and do not replace user training. Intensive instruction and practice is essential for the use and maintenance of the rescue harness with carrying function. In addition, users must be physically and mentally fit and not under the influence of intoxicants. Inadequate training, incorrect use and misuse of the product can lead to accidents. The instructions for use and the information and instructions contained therein must be carefully observed and complied with.

1. General application

The rescue harness with carrying function has been developed for rescuing injured persons from areas with difficult access. The rescue harness is used for the horizontal and vertical transport of injured persons from exposed positions for rope and air rescue. The injured person is secured with the integrated rescue harness and additionally wrapped in a flexible plastic plate using fixation strap fasteners. The aim is to ensure the safety of the injured person and minimise the physical exertion of the rescue team. The rescue harness with carrying function may only be used within the specified operating conditions and for the intended purpose.

Any combination of different equipment can cause unforeseeable hazardous situations and have a negative impact on the safety of the user.

2. Areas of application

Rescue from confined spaces, rescue from heights and depths, rope and air rescue, rescue in hard-to-reach areas, industrial accidents, evacuation from natural disasters, emergencies in urban environments, military or tactical operations and road accidents

3. User training

The user must be sufficiently trained and familiar with the use of the rescue stretcher before using it for the first time. In addition, repeated instruction in these instructions for use is necessary during the introduction and subsequent further training.


The training course must be documented in detail and repeated at regular intervals in accordance with applicable standards and national regulations. The type, scope, and date of the training must be recorded chronologically.

4. Product description


The RsqTec TURTLE FrameX is a versatile and certified rescue system that offers a combination of patient transport equipment and rescue harness. It has been tested in accordance with the standards DIN EN 1497:2007 and DIN EN 1865-1:2015 clause 5.4.1 and is particularly suitable for the rescue and transport of injured or accident victims. The application starts with opening all buckles and placing the straps to the sides to allow free handling. The adaptive leg support is inserted into the square tube of the rescue stretcher and locked in the integrated locking system. The leg support can also be retrofitted during transport. The injured person is carefully positioned on the rescue stretcher, with the top of their head flush with the upper edge of the protective tray. The leg straps are then correctly attached by passing the belt under the thigh and wrapping it around the leg at the crotch. The buckles are audibly and perceptibly engaged and the straps are pulled tight, but not constricting, so that the person's position remains stable. The shoulder strap is passed over the head and placed on the chest. Here too, the buckles are connected and the straps are tightened slightly to ensure that the chest remains free to move for breathing. The head fixation made of non-slip fabric is attached to the person's forehead or helmet and secured with adjustable side straps. The lashing straps are connected crosswise and pulled taut, to securely shape the protective tray of the stretcher. A head hanger is used for suspension, which is passed over the person and

fastened to the stretcher's rescue harness with side buckles. The black elasticated loops and the blue straps must be correctly aligned. The person's arms are attached to the elasticated strap loops and the legs are secured to the leg support with Velcro fasteners. In addition, an adjustable connecting sling can be connected to the karabiner of the lifting device to adjust the transport position from horizontal to vertical and vice versa. The rescue stretcher is now ready for transport. The O-ring, as the main attachment point, is used to attach and transport the stretcher. Other optional attachment points are labelled with an "A" and include a D-ring in the chest strap and an attachment loop in the centre of the head harness, the use of which requires the disconnection of the buckles from step 8. Thanks to these well-thought-out functions and the optional adjustment options, the rescue stretcher ensures safe and flexible transport of injured persons, allowing both horizontal and vertical positions.


5. Photographic representation

6. Step by step instructions - RsqTec TURTLE FrameX

Step	Picture	Description of the
1		<p>Open the plug-in fasteners: Fully open all plug-in fasteners on the appliance.</p> <p>Position the straps: Place all straps to the sides so that they are not in the way.</p>
2		<p>Insert the leg support: Insert the adaptive leg support into the square tube of the rescue stretcher.</p> <p>Locking: Lock the leg support with the integrated locking system.</p> <p>Subsequent attachment: Leg support can also be attached during transport.</p>
3		<p>Positioning Person: Carefully lay the injured person on their back on the rescue stretcher. Ensure that the top of the person's head is flush with the upper edge of the protective tray.</p>
4		<p>Pull the belt through: Guide the long side of the belt under the respective thigh from the outside to the inside.</p> <p>Lay the loops: Place around the leg at the crotch.</p> <p>Connect: Allow the plug-in fasteners to engage audibly and tactilely. Position the padding pads on the push-fit fasteners and between the legs</p> <p>Pull tight: Tighten the straps (not constricting) and maintain the person's position.</p>

<p>5</p>		<p>Position the shoulder strap: Guide over the head and place on the chest.</p> <p>Connect the buckles: Connect the shoulder strap in front of the chest to the side chest straps.</p> <p>Tightening: Tighten the straps slightly (do not constrict).</p> <p>Check breathing: Ensure freedom of movement of the chest for breathing.</p> <p>Check the head position: Check the head position from step 2.</p>
<p>6</p>		<p>Position the head fixation: Position the non-slip fabric head fixation on the forehead/helmet.</p> <p>Fasten the side straps: Fasten the adjustable straps at the side.</p> <p>Medium adjustment: Adjust the flexible black strap in the centre behind the head.</p>
<p>7</p>		<p>Apply lashing straps: Connect the fasteners crosswise.</p> <p>Tighten: Tighten the straps (not constricting), bring the protective tray into shape</p>
<p>8</p>		<p>Position the suspension: Guide over the person with the headband.</p> <p>Connect the buckles: Connect the side buckles to the rescue harness of the stretcher.</p> <p>Check alignment: Black elasticated loops on the outside, blue straps must not be twisted.</p>
<p>9</p>		<p>Use the elasticated strap loops: Pass your hands through the loops.</p> <p>Hook in the arms: Insert the arms into the elasticated strap loops.</p> <p>Check position: Align the loops on the blue side straps towards the outside.</p>

<p>10</p>		<p>Fix the legs: Attach the legs to the leg rest using the Velcro fastener.</p> <p>Optional connection: Connect the adjustable connecting sling to the karabiner of the lifting device.</p> <p>Adjustability: Sling makes it easier to adjust the transport position from horizontal to vertical and vice versa.</p>
<p>11</p>		<p>Fix the legs: Attach the legs to the leg rest using the Velcro fastener.</p> <p>Preparing for transport: The rescue stretcher is now prepared for transport.</p> <p>O-ring (main attachment-point): The O-ring in the illustration is used for attachment and transport.</p> <p>Optional attachment points: These are labelled with an "A". 1st D-ring, located in the chest strap (see step 5). 2nd stop loop is located in the centre of the headband. To use this loop, the fasteners from step 8 must be separated.</p> <p>Optional sling: The adjustable connecting sling can (but does not have to) be connected to the carabiner of the lifting device in order to change the transport position.</p>

7. Attachment and transporting - Applications

The following application options may only be used after the harness straps of the rescue stretcher have been fully fastened (see chapter 6).

Requirements:

- The anchor point must always be above the rescue stretcher.
- A minimum load-bearing capacity of the anchor point of 12 kN is required.
- All attachment points on the rescue stretcher are labelled with an "A" symbol.

Possible applications:

1. O-ring (main attachment point)
 - Horizontal transport: The snap carabiner is connected to the O-ring.
 - Vertical transport: The snap carabiner is disconnected from the O-ring.

Advantage: A change from the horizontal to the vertical position can be carried out under load. This is done by relieving the load on the rescue stretcher using the leg support, which allows the snap carabiner to be released from the O-ring.



2. Loop on the headband (optional)
 - Application: Vertical transport.
3. D-ring in the rescue harness (optional)
 - For additional attachment options.

Note: Always observe the safety instructions and the maximum permissible loads during use.

8. Technical data and parameters

Category	details
Application	RsqTec Turtle FrameX
Dimensions	In use: 1200 x 650 mm In use with leg rest: 1870 x 650 mm
Nominal load	140 kg
Weight (net)	12.5 kg
Certification	DIN EN 1497:2007 (Rescue harness) DIN EN 1865-1:2015 (patient transport equipment for rescue vehicles) ANSI/ASSPZ359.12-2019 and CSA Z259.12-16 (quick-release buckles)
Lifetime	10 years
Temperature range	-40°C to +60°C
Areas of application	Rescue from confined spaces Rescue from heights and depths Rescue by rope or from the air Rescue in areas that are difficult to access Industrial accidents Evacuation in the event of natural disasters Emergencies in urban environments Military or tactical operations Traffic accidents
Special features	Mounting options: horizontal, vertical and in the rescue harness Protection: Equipped with a protective shell Handling: Four side carrying handles for easy handling Stability: Integrated reinforcement frame ("FrameX")

9. Description of the label - nameplate

label	description	nameplate
Mittelmann Sicherheitstechnik	Manufacturer: Mittelmann Sicherheitstechnik GmbH & Co. KG	 <p>Mittelmann Sicherheitstechnik Bessemerstr. 25, 42551 Velbert www.mittelmann.com RsqTec TURTLE FrameX DIN EN 1497:2007 DIN EN 1865-1:2015 Ziffer 5.4.1 Nennlast: 140 kg rated load Made in Germany CE 0158 Baujahr: xx/xxxx year of manufacture   Fabr.-Nr.: xxxxxxxx serial number <input type="text"/></p>
Bessemerstrasse 25, 42551 Velbert	Postal address of the manufacturer	
www.mittelmann.com	Internet address of the manufacturer	
RsqTec	Product designation	
TURTLE FrameX	Type designation	
DIN EN 1497:2007 DIN EN 1865-1:2015 Clause 5.4.1	Numbers and years of the standard documents to which the equipment conforms	
Rated load: 140 kg	Specification of the maximum permissible body weight	
Made in Germany	Designation of origin	
CE 0158	Identification number of the quality monitoring notified inspection body	
Year of construction: xx/xxxx	Month / year of manufacture	
	Symbol to indicate that the instructions for use must be observed	
Serial no.: xxxxxxxxxx	Serial number	

10. Material specifications

The product consists of the following materials:

- Polyamide and polyester webbing
- Polyamide fabric cuttings
- Polyester sewing thread
- Fittings made of aluminum, steel and stainless steel
- Plastic parts
- Aluminum and stainless-steel push-in fasteners
- Aluminum or stainless-steel snap hooks

11. Storage and transport

- Any kind of contact with acids, corrosive liquids and oils must be avoided at the storage location.
- If the product has become damp or wet, it must be dried naturally before storage. Always dry naturally, never in the vicinity of fire or other sources of heat.
- The product should be placed in a suitable transport bag for transport.

12. Operational test

The inspection points below are to be checked before each use and apply only to stretchers that are not vacuum-sealed:

- Webbing or fabric damage due to abrasion, cutting, heat or chemicals
- Any type of damage to seam connections
- Any type of damage to the supporting textile structure
- Any type of soiling caused by technical liquids or solids that can no longer be cleaned
- Damage to the Velcro fasteners
- Deformation or corrosion of the fittings or push-in fasteners
- Defective fasteners or buckles
- Legibility of product labelling

As soon as one of the characteristics is determined to be "not OK" during the pre-use inspection, the rescue harness with carrying function must not be used.

The inspection points are not required for sealed stretchers, provided that the annual inspection of the vacuum unit has been inspected each year.

13. Safety requirements

- Any combination of different equipment can cause unforeseen hazardous situations and negatively affect the safety of the user
- All instructions in this user manual must be observed to protect the user and the rescue harness with carrying function. The rating plate must be legible.
- To avoid the occurrence of suspension trauma to the injured person, the user must have an emergency plan in place. (Suspension trauma can occur if the return flow of blood from the legs is impeded or lost during prolonged, motionless hanging in a harness. This can lead to orthostatic circulatory shock. The injured person cannot perform any muscle pumping or other relieving movements due to the fixation by the textile elements. The duration of hanging in a vertical position must therefore also be minimized.
- It is necessary to check the attachment to the person again before use.

- For the safety of the person to be rescued, it is essential that the necessary clearance for the rescue procedure is ensured before each use of the rescue harness with carrying function, so that no impact with the ground or other obstacles is possible.
- Hazards such as extreme temperatures, exposure to sharp edges, cuts, abrasion, slack ropes, chemicals, electrical influences or climatic effects can impair the function of the equipment

14. Inspection and general information

- If there is any doubt about the safe condition of the system or its components, they must be replaced immediately. This must be carried out by the manufacturer or another competent person.
- Systems damaged by a fall must be withdrawn from use and returned to the manufacturer or an authorized repair workshop for maintenance or repair and retesting.
- **Do not make any changes to the product.**
- **Repairs may only be carried out by the manufacturer.**   
- An instruction manual / inspection book is supplied with every rescue harness with carrying function. This must be kept with the equipment 
- The contents of the instructions for use of the products used in conjunction with the rescue harness with carrying function and the applicable accident prevention regulations must be observed.
- The declaration of conformity can be downloaded from the Internet at www.mittelmann.com.
- As required, but at least once a year, an inspection must be carried out by the manufacturer or by an expert (see point 12 Operational test). This inspection must be documented and dated in the inspection log provided. In addition, the next inspection date must be visibly noted on the type plate (e.g. with an inspection sticker). 
- The rating plate with date of manufacture and serial number is located in the rating plate pocket with the i symbol on the left shoulder strap.
- If the equipment is resold to another country, the reseller must provide instructions for use, maintenance, periodic inspection and repair in the language of the other country.
- The rescue harness with carrying function must not be used for fall arrest purposes.
- The rescue harness with carrying function may only be used and/or stored unprotected in the offshore area for a limited period of time. The metal components must be checked for signs of corrosion with appropriate care and at shorter intervals.

15. Cleaning instructions

- After each use, the rescue harness with carrying function must be thoroughly cleaned and then stored dry in a ventilated room protected from UV radiation.
- Avoid contact with acids, corrosive liquids and oils. If unavoidable, clean as quickly as possible with mild detergent and hand wash in plenty of water (approx. 30 °C). Always dry naturally, never near fire or other sources of heat.
- Do not dry clean, tumble dry or spin dry and do not use bleach.  
- If you have any questions about disinfection, please contact the manufacturer.
- Any body fluids on the rescue harness with carrying function may only be cleaned with water.

16. Service life

The service life of the TURTLE FrameX rescue stretcher depends on the respective operating conditions. Based on previous experience, a service life of 10 years can be expected under normal operating conditions. The year of manufacture of the rescue stretcher is indicated on the nameplate. If stored correctly and the sealed packaging is checked annually in accordance with the specified guidelines, the rescue stretcher can be stored in the vacuum-sealed packaging for up to 15 years, with a maximum service life of 15 years. The rescue harness with carrying function must be returned to the manufacturer in the event of any type of damage to any components for the damage to be analyzed and, if necessary, repaired. Continued use of the rescue harness with damaged components will result in immediate loss of the liability claim.

Notified body involved in production monitoring:

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