Mittelmann

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
RsqTec **TURTLE**

DIN EN 1497:2007

DIN EN 1865-1:2015 Ziffer 5.4.1

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DE

Instructions for use / test book

GB

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Instructions for Use

These instructions for use are general instructions for the product described and are not a substitute for user training. Intensive instruction and practice in the use and maintenance of the rescue harness with carrying function is mandatory, in addition, the user must also be both physically and mentally fit. Inadequate training, incorrect use or misuse of the product can lead to accidents. The advice and instructions contained in the instructions for use must be carefully taken note of and complied with.

Foreword

The strengths of Mittelmann Sicherheitstechnik and of Deutsche Windtechnik have been brought together for the development of the "TURTLE" rescue harness with carrying function. Hence, it was possible to combine development know-how for special rescue equipment together with the industry-specific requirements and ideas of the wind energy sector.

During the development time, the rescue harness with carrying function has been tested in regular training operations on wind turbines of 5 different manufacturers and improved continuously throughout this phase.

1. General Application

The rescue harness with carrying function has been developed for the rescue of injured persons from difficult to access areas. The rescue system enables horizontal and vertical transportation of injured people from exposed locations for rope and air rescue. The encasing of the person takes place by means of tear-resistant and water-repellent nylon fabric and the injured person is secured with an integrated rescue harness. The objective is to ensure the safety of the person and to minimize the physical exertion of involved rescue workers as much as possible. The rescue harness with carrying function must only be used within the specified operating conditions and for its intended purpose.

Through interaction, every combination of different equipment can give rise to unforeseen hazardous situations and negatively affect the safety of the user.

2. Application Areas

Rescue from wind turbines, confined spaces and rescue from heights.

3. User Training

The emergency personnel must be adequately trained and familiar with the use of the rescue harness with carrying function before its first use. In addition, in the introduction and the subsequent further training, the repeated teaching of these instructions for use is necessary.

The training course must be strictly documented and be repeated at least once per year. The type, the scope and the date of the training must be chronologically recorded.

4. Product Description

The “*TURTLE*” rescue harness with carrying function is a combination of patient transport means and a rescue harness and has been tested and certified in accordance with DIN EN 1497:2007 and DIN EN 1865-1:2015 Ziffer 5.4.1.

TURTLE is a rescue system, made using textile materials, for the rescue and transport of accident victims. The textile sheath (a) encases the person from the chest to the buttocks area. Fastening is done by means of large size Velcro strips. Two lashing belts (b) with adjustable Cobra snap buckles serve for tightening of the protective cover around the person's body. A system of supporting longitudinal, lateral and diagonal belts inside the stretcher carry the body weight of the person during the rescue and transport operations. The head strap (c) on the upper part of the stretcher serves for the pulling it horizontally. Opposite to this there is a buttocks strap (d). The shoulder straps (e) are adaptable by means of adjustable buckles that in turn are connected to the horizontal chest strap (f) via two snap buckles. Additional supporting straps are the leg straps left (g) and right (h) that also can be adjusted in length for each person by means of the Cobra snap buckles. The rescue and transport of a horizontally lying person is performed via anchoring to the four suspension loops (i) that can be optionally equipped with two round rings (two loops per ring) or with four shackles. The transport of a vertically hanging person takes place by means of the vertical suspension loop (k). In addition, there are carrying loops with handles (j) on the outside of the rescue harness for manual transport.

5. Designation of the Individual Parts

No.	Designation	Anchor points for these lying positions
(a)	Textile sheath (blue)	
(b)	2x lashing belts (black)	
(c)	1x head strap (yellow + blue belt)	For horizontal pulling
(d)	1x buttocks strap (yellow + blue belt)	For horizontal pulling
(e)	2x shoulder straps (blue)	
(f)	1x chest strap (blue)	
(g)	1x leg strap left (red belt)	
(h)	1x leg strap right (green belt)	
(i)	4x suspension loops (belt straps with round rings / shackles)	For horizontal lifting
(j)	4x carrying loops with handles (black + blue)	No anchor point
(k)	1x suspension loop	For vertical lifting

5.1 Illustrative Representation



Illustration Julia Kotulla

5.2 Photographic Illustration

6. Belt Fitting of the RsqTec “TURTLE”

Before starting the belt fitting, all belts and snap buckles are opened and laid to one side. Now the injured person is laid on his or her back on the RsqTec *TURTLE* so that his or her crown is approx. 5 cm away below (caudal) from the upper (cranial) edge of the protective tray. Then the tabularly represented step-by-step instructions for the correct belt fitting are carried out.

Illustration	Description
	Coming from the side, pass the long belt side of the leg strap (red (g) and green (h)) under the respective thigh of the person. Guide them up through the crotch and connect them to the relevant snap buckles (the snap buckles must audibly and tactilely snap into place). The strap is pulled tight (not constricting) without changing the position of the person in the process.
	Lay the shoulder strap (e) over the head of the person in the longitudinal direction, guide it over the shoulder and connect it at the front of the chest with the side parts of the chest strap (f). The strap is pulled tight (not constricting). The complete freedom of movement of the rib cage for breathing must be strictly maintained.

Before the enclosure of the upper body, **ensure** that the head position on the head strap (c) and the strapping are comparable with that shown in the illustration.


The head fixation consists of an adhesive textile fabric and is fixed to the forehead of the injured person or to the front surface of the helmet with the red and green webbing. The position of the head restraint can be adjusted using the flexible black strap.

Lay the textile sheath closely around the person and close it using the Velcro fastener.

	Connect the protective tray with the aid of the lashing belts and snap buckles and pull them close around the body of the injured person (<u>not constricting</u>).
	The person is now prepared for pulling and lifting with the <i>rescue harness TURTLE</i> in various positions

Accessories	Description
	Foot strap for load relief of the legs in vertical hanging positions. Lowers the risk of orthostatic shock. The foot sling or foot bag is attached to the black webbing loops below the leg loops by means of push-in fasteners.

	The backpack for the RsqTec TURTLE can be used as a footmuff (only for horizontal transport) with the help of the stiffening plate supplied. To do this, insert the stiffening plate with the narrow plastic side into the viewing window. The slight bend of the plate points upwards. The footmuff is then attached to the rescue stretcher by means of the integrated buckles (see attachment of foot sling). The footmuff is also connected to the round rings with an adjustable lanyard (not included).
	The neck support pillow is intended to relieve tension and is placed on the red fleece band in the head area of the rescue stretcher and thus connected.

7. Pulling and Lifting in Various Positions

The following possible applications first take place after the successful carrying out of the belt fitting of the RsqTec "TURTLE" (see Chapter 6.)

The anchor point for all positions must be lie above the user and guarantee a minimum load carrying capacity of 10 kN.

Horizontal lifting of the person with the RsqTec TURTLE:

- (i) 4x suspension loops (belt straps with round rings or shackles)

The four suspension loops (i) are removed from their side packing bag and attached to the lifting tackle of the hoist.

Lift the person with the aid of a suitable hoist, whereby the horizontal position of the person must be constantly checked and, if necessary, be corrected by changing the lying position of the person.

Vertical lifting of the person with the RsqTec *TURTLE*:

Note: In order to prevent injury to the person, the person must only be lifted in the vertical position for a short time (for a few minutes).

- (k) 1x suspension loop

Attach the lifting tackle of the hoist to the suspension loop (k) (centred, within the loop provided for this) and lift the person with the aid of a suitable hoist.

During this process, special attention must be paid so that no excessive pressures occur for the person. Where necessary, the position of the relevant belt must be checked.

Horizontal pulling of the person with the RsqTec *TURTLE*:

- (c) Head strap

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With the aid of the head strap (c), the person can be pulled in the direction of the head horizontally and up to an inclination of approx. 45°. The position of the pulling must take place approximately in the plane of the longitudinal axis of the stretcher and can be done by hand or with the aid of suitable lifting equipment.

The same process can also take place in the direction of the feet of the person with the aid of the buttocks strap (d).

8. Technical Data and Parameters

Physical Dimensions

Maximum length	1050 mm = 1.05 m
Maximum width	650 mm = 0.65 m
Weight without a person	approx. 7.5 kg

Load Limits


Maximum nominal load	140 kg
----------------------	--------

Temperature Limits

Minimum outdoor temperature	- 40 °C = - 40 °F
Maximum outdoor temperature	+ 60 °C = + 140 °F

9. Meaning of the Marking

Mittelmann Sicherheitstechnik	Manufacturer Mittelmann Sicherheitstechnik GmbH & Co. KG
Bessemerstraße 25, D-42551 Velbert	Postal address of the manufacturer
www.mittelmann.com	Internet address of the manufacturer
RsqTec	Product designation
TURTLE	Type designation
DIN EN 1497:2007 DIN EN 1865-1:2015 Ziffer 5.4.1	Number and year of the document with which the equipment complies
Nominal load: 140 kg	Specification of the maximum permissible body weight
Made in Germany	Designation of origin
 0158	Identification number of the quality monitoring notified inspection body
Year of manufacture: xx/xxxx	Month / year of the manufacture
	Symbol for indicating that the instructions for use must be observed
Serial No.: xxxxxxxxxx	Serial number

10. Material Specifications

The RsqTec *TURTLE* consists of the following materials:

- Straps in polyamide
- Material sections in polyamide
- Sewing thread in polyester
- Fitting components in aluminium, steel and stainless steel
- Plastic components in polyethylenes
- Snap buckles in aluminium and stainless steel
- Plastic parts in synthetic material
- Carabiner hooks in aluminium or stainless steel

11. Storage and Transport

- The RsqTec *TURTLE* must be stored dry, protected from direct sunlight and at room temperature (approx. 17 °C – 23 °C).

- At the storage location, all types of contact with acids, corrosive liquids and oils must be prevented.

- Whenever the RsqTec *TURTLE* becomes wet, it must be dried naturally before being stored.
- During the transport of the RsqTec *TURTLE*, it must be stored in a suitable transport bag.

12. Damage Analysis

Before each use of the rescue harness with carrying function it must be checked for the following signs of wear and / or damage:

- Strap or textile damage due to abrasion, cuts, heat or chemicals.
- Damage to any type of seam
- All types of damage to the load-bearing textile structure.
- All types of contamination through technical liquids that can no longer be cleaned.
- Damage to the Velcro fasteners
- Deformation or corrosion of the fitting components or snap buckles
- Defective snap buckles or clasps
- Check the legibility of the product marking.

Immediately one of the damage analysis points described above is detected before use, the rescue harness with carrying function must not be used.

13. Safety Requirements

- Through interaction, every combination of different equipment can give rise to unforeseen hazardous situations and negatively affect the safety of the user.
- All instructions for the protection of the user and the rescue harness with carrying function given in these instructions for use must be strictly observed. The rating plate must always be visible and legible.
- In order to avoid a suspension trauma, an emergency plan on the part of the user must be available for every type of use of the rescue harness with carrying function. A suspension trauma is a life-threatening state of shock that occurs after long, motionless hanging in a harness system. This can lead to unconsciousness and to cardiac arrest.
- This emergency plan must ensure that the vertical hanging time of the person in the rescue harness is kept as short as possible and so prevent massive damage to health.
- It is vital that the belt fitting to the person is regularly checked during its use.

- For the safety of the person to be rescued, it is essential that before every use of the rescue harness with carrying function the free space necessary for the rescue operation is ensured so that no impact on the ground or with other obstacles is possible.
- Hazards, such as for example, extreme temperature, stress from sharp edges, cuts, abrasion, slack rope formation, chemical effects, electrical influences or climatic effects can adversely affect the function of the equipment.

14. Inspection and General Instructions

- Should there be doubts about the safe condition of the system or the component parts, these must be replaced immediately. This must be carried out by the manufacturer or another qualified person.
- Systems stressed by a fall must be taken out of use and be returned to the manufacturer or a qualified workshop for servicing and re-examination.


- **Do not make any changes to the product.** 
- **Repairs must only be carried out by the manufacturer.**  
- The instructions for use / test book are delivered with every rescue harness with carrying function. These must be kept with the equipment. 
- Strictly observe the instructions for use of the products used in connection with the rescue harness with carrying function as well as the relevant applicable accident prevention regulations.
- The Declaration of Conformity can be downloaded from the www.mittelmann.com page.
- As required, but at least once per year, an inspection by the manufacturer or a qualified person must be carried out (see also Point 12 Damage Analysis). This inspection with the date must be documented in the test book supplied. In addition, the next inspection date must be clearly recorded on the rating plate.
- The rating plate with the year of manufacture and the serial number can be found in the rating plate pocket on the left shoulder strap (e).
- If the equipment is sold on in another country, the re-seller must provide instructions for the use, the maintenance, the regular inspection and repair in the language of the other country.
- The rescue harness with carrying function must not be used for fall arrest purposes.
- The rescue harness with carrying function must only be used only for a shortened period of time in offshore environments. The metal components must be inspected for corrosion with increased carefulness and in shorter time intervals.

15. Cleaning Instructions

- The rescue harness with carrying function must be thoroughly cleaned after every use and be stored dry in an airy and shady space.

- Avoid contact with acids, corrosive liquids and oils. If unavoidable, clean as soon as possible by hand washing with mild detergent and plenty of water (approx. 30°C). Always allow to dry naturally, never in the vicinity of fire or other sources of heat.

- Please do not use dry cleaning, tumble dryers, spin dryers or bleaching processes.

- For questions regarding disinfection, please contact the manufacturer.
- Any body fluids must be cleaned off using water only.

16. Service Life

The service life is dependent on the respective usage conditions. On the basis of past experience, under normal usage conditions the service life of the *TURTLE* rescue harness with carrying function can be assumed to be 10 years. The year of manufacture of the rescue harness is recorded on the rating plate.

In the event of any kind of damage to any component, the rescue harness with carrying function must be available to the manufacturer for the evaluation of the damage. The continued use of the rescue harness with damaged components shall lead to the immediate loss of any and all rights to liability claims.

Notified body involved in the production monitoring:

DEKRA Testing and Certification GmbH

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